



# 陕西绿清生物工程有限公司

## Shaanxi Green Bio-Engineering Co.,Ltd

### Determination of Luteolin Purity

A. 1 Methods: The samples are extracted by ultrasonic wave and detected by high performance liquid chromatography (HPLC).

A.2 Test instruments and appliances

A.2. 1 Analytical balance with accuracy of 0.00001g.

A.2.2 Ultrasonic cleaner: 250W, 20kHz

A.2.3 High performance liquid chromatograph

A.3 Reagents and solutions

A.3. 1 Acetonitrile, analytical pure water

A.3.2 Secondary distilled water

A.3.3 Luteolin control product.

A.3.4 Preparation of mobile phase: mixed with acetonitrile-water (30:70 ratio) and obtained by filtration with microporous membrane.

A.3.5 Detector and detection wavelength: UV spectrophotometer with wavelength of detection 360nm.

A.4 Methods of operation

A.4. 1 Preparation of control solution: Luteolin control product (accurate to 0.01 mg) is accurately weighed and added with methanol to prepare a solution containing 70 $\mu$ g per 1mL as control product solution.

A.4.2 Preparation of test solution: Take A sample of luteolin (about 10mg), weigh it accurately, dissolve it with ultrasonic methanol, and use it as the test solution.

A.4.3 Determination methods

Precisely absorb 10 $\mu$ l of control solution and test solution respectively, inject them into the liquid chromatograph for determination.

#### A.5 Calculation of results

The purity of luteolin is calculated according to Equation (B. 1) :

$$\text{Luteolin (\%)} = \frac{S_1 \times C \times A}{S_0 \times (M - M \times B)} \times 100\% \dots\dots\dots (\text{B.1})$$

Where,

S<sub>1</sub>-- Peak area value of the test product solution;

S<sub>0</sub>-- Peak area value of the control product solution;

C-- Concentration of the control product solution (mg/mL);

A-- Control product purity (%);

B-- Moisture purity of the test product (%);

M-- Concentration of test product solution (mg/mL).